



Abutment Screw Torque and Driver

Manufacturer	Torque N-cm							Driver	
	15	20	24	25	30	35	45	Titanium Screw	Gold Screw
Atlantis™ titanium and zirconia abutments utilize the same torque and driver setting as the original implant manufacturer									
Astra 3.0	X							0.050"- 1.25mm hex	
Astra 3.5-4.0		X						0.050"- 1.25mm hex	
Astra 4.5-5.0				X				0.050"- 1.25mm hex	
BioHorizons® External, Internal, Internal Tapered					X			0.050"- 1.25mm hex	
Biomet 3i™ External Hex						X		0.048" hex	Square
Biomet 3i™ Certain® *		X						0.048" hex	
BlueSkyBio One Stage					X			Star & 0.048" hex	
BlueSkyBio Trilobe					X			Unigrip & 0.048" hex	
BlueSkyBio Internal Hex					X			0.048" hex	
Brånemark (Nobel Biocare)						X		Unigrip	Square
Camlog™ & Conelog®		X						0.050"- 1.25mm hex	
Dentium Super Line & Implantium					X			1.28mm hex	
Dentsply Friadent			X					0.9mm hex	
Dentsply Ankylos® C/X	X							1mm hex	
Glidewell iNCLUSIVE® Tapered 3.7, 4.7, 5.2						X		0.050"- 1.25mm hex	
Hiossen/Osstem HG mini (3.5)		X						1.2mm hex	
Hiossen/Osstem HG standard (4.0, 4.5, 5.0)					X			1.2mm hex	
Imtec® 3M Endure™						X		0.050"- 1.25mm hex	
Implant Direct™					X			0.050"-1.25mm hex	
Keystone/Lifecore Genesis* & Prima™ *					X			Square	
Keystone/Lifecore Renova®					X			0.048" hex	
Keystone/Lifecore Restore®					X			0.048" hex	Square
Mega'Gen EZ Plus						X		0.048" hex	
Mega'Gen Rescue							X	0.048" hex	
Mega'Gen ExFeel internal (3.5)				X				1.2mm hex	
Mega'Gen ExFeel internal (4.1, 4.8, 5.5)						X		1.2mm hex	
Mega'Gen ExFeel external (3.3)				X				1.2mm hex	
Mega'Gen ExFeel external (3.75, 4.0,4.5, 5.0, 5.5)						X		1.2mm hex	
MIS® Biocom & Seven					X			0.050"-1.25mm hex	